

PCD-41D

Pockels Cell Driver with two-step pulse and voltage doubling


Generation of a two-step high voltage pulse with nanosecond rising edges.

First step duration of the output high voltage pulse is determined by trigger pulse duration.


Second step amplitude of the output HV pulse is equal to the double voltage of first step.

Driver can be used for ultrafast optical beam modulation and deflection,

control of laser regenerative amplifiers and pulse pickers.



- Two-step high voltage pulse;
- Nanosecond rising edges;
- Controlled pulse duration;
- Pulse jitter is less than 200 ps;
- Built-in pulsed high voltage source;
- Additional output signal for HV pulse amplitude measurement;
- Pulse amplitude adjustment by trimmer or external analog signal;
- Low voltage power supply 12 V.



Specifications


Output amplitude of the first step ¹	2800 ÷ 3600 V
Output amplitude of the second step (double voltage of first step) ¹	5600 ÷ 7200 V
Unevenness of the amplitude on the plateau	< 5 %
Output voltage pulse-to-pulse instability	< 1 %
Max pulsed current	20 A
Optimal load capacitance	4 ÷ 12 pF
HV pulse rise time of first step ²	2 ± 1 ns
HV pulse rise time of second step ²	2 ± 1 ns
Pulse duration of the first step ³	2 ÷ 600 ns
Duration of the second step plateau	≈ 1000 ns
Max HV pulse repetition rate ²	1000 Hz
Trigger voltage (input impedance is 50 Ω)	5 ± 0.5 V
Output pulse delay vs. trigger pulse ⁴	16 ± 3 ns
HV pulse jitter	< 0.2 ns
External power supply voltage	12 ± 4 VDC
External power supply current	500 mA
Operating temperature range	-40 ÷ +60 °C
Dimensions	80 × 50 × 15 mm
Mounting hole pattern (Ø 3.2 mm)	74 × 44 mm
Weight (OEM version)	100 g

- 1 HV amplitude is controlled by built-in trimmer RV1. The range of voltage corresponds to a specific driver model.
- 2 Depends on the load impedance and output voltage amplitude.
- 3 Is determined by an external trigger-generator.
- 4 Delay depends on the trigger pulse and HV output amplitudes. The higher the HV output amplitude, the shorter delay.

Leading the Light 1



Connection diagram



Connector X1 (input) - 53261-0471 (Molex):

Pin 1 (red) Power supply + 12 V & 500 mA max;

Pin 2 (black) Power supply GND;

3 Pin 3 (green) Output voltage measure signal;

DC voltage scale 1:10000; Note 1

4 Pin 4 (black) Output voltage measure GND.

Connector J2 (output) - SM02B-BHSS (JST):

1 Pin 1 (HV red) First high voltage output;

2 Pin 2 (HV red) Second high voltage output.

Connector X3, X4 (input) - SMA-J-P-H-ST-EM1 (Samtec):

Trigger pulse input +5 V; impedance 50 Ω .

Connector X3 input trigger pulse determines output HV pulse duration for first step.


Connector X4 not available.

Note 1

Pin 3 and **Pin 4** can be can be used for setting the amplitude of output voltage pulse from –20 % to +2 % (relative to value have been set by trimmer RV1).


If 0 V is set on Pin 3 from from an external source, the pulse amplitude will be ~2 % higher than the set value.

If **5 V** is set on the **Pin 3**, the pulse amplitude will be lower by ≈ 20 %. The input impedance of the **Pin 3** is $45 \text{ k}\Omega$.




2


Leading the Light



Output waveform of PCD-41D - two-step pulse with voltage doubling





Block diagram of PCD-41D - two-step pulse with voltage doubling



Waveforms of the transmitted light beam

Pulse voltage: $V_{OUT} = 2 \times 3600 \text{ V}$; Pockels cell half-wave voltage: $V_{\chi/2} = 3800 \text{ V}$





4

Leading the Light



Leading the Light

www.teo.technology info@teo.technology

in LinkedIn

